


Raman spectroscopy

Raman spectroscopy, being a spectroscopy, is a technique that uses the interaction of light with matter to analyze the chemical composition and structure of materials. **Raman spectroscopy is based on the principle that when a beam of light (usually a laser) is scattered by a molecule, some of the photons undergo a change in their energy and wavelength.** This change, called the Raman shift, is proportional to the vibrational frequency of the molecule and can be used to identify its molecular fingerprint. Raman spectroscopy can provide information about the molecular bonds, functional groups, crystal structure, phase transitions, and stress/strain of materials. Raman spectroscopy is a non-destructive and versatile technique that can be applied to solids, liquids, gases, and biological samples.

Raman Spectroscopy


1 in 10^7 photons is scattered inelastically



Raman spectroscopy. Principles: Light scattering

Light scattering is a physical phenomenon that occurs when **electromagnetic waves interact with particles or molecules in a medium. The scattering process can change the direction, frequency, polarization, phase, or intensity of the incident light.**



Depending on the size and shape of the scatterers, as well as the wavelength of the light, different types of scattering can occur, such as Rayleigh scattering, Mie scattering, Tyndall scattering, and Raman scattering. Light scattering is important for many applications in science and engineering, such as atmospheric optics, spectroscopy, microscopy, optical communication, and remote sensing.



Raman spectroscopy. Principles

When light is scattered by molecule, the oscillating electromagnetic field of a photon induces a polarisation of the molecular electron cloud **which leaves the molecule in a higher energy state with the energy of the photon transferred to the molecule**. This can be considered as the formation of a very short-lived complex between the photon and molecule which is commonly called the **virtual state of the molecule**. The virtual state is not stable, and the photon is re-emitted almost immediately, **as scattered light**.


In the vast majority of scattering events, the energy of the molecule is unchanged after its interaction with the photon; and the energy, and therefore the wavelength, of the scattered photon is equal to that of the incident photon. This is called **elastic (energy of scattering particle is conserved) or Rayleigh scattering** and is the dominant process.



Raman spectroscopy. Principles


In a much rarer event (approximately 1 in 10 million photons)² Raman scattering occurs, which is an **inelastic scattering process with a transfer of energy between the molecule and scattered photon**. If the molecule **gains** energy from the photon during the scattering (excited to a higher vibrational level) **then the scattered photon loses energy and its wavelength increases which is called Stokes Raman scattering**

Inversely, if the molecule **loses** energy by relaxing to a lower vibrational level **the scattered photon gains the corresponding energy and its wavelength decreases; which is called Anti-Stokes Raman scattering**. Quantum mechanically Stokes and Anti-Stokes are equally likely processes. However, with an ensemble of molecules, the majority of molecules will be in the ground vibrational level (Boltzmann distribution) and Stokes scatter is the statistically more probable process. As a result, **the Stokes Raman scatter is always more intense than the anti-Stokes and for this reason, it is nearly always the Stokes Raman scatter that is measured in Raman spectroscopy.**



Raman Spectrum of CCl_4

CCl_4 is a tetrahedral molecule with three pronounced Raman active vibrational modes in the 100 cm^{-1} to 500 cm^{-1} wavenumber region (there is an additional peak at $\sim 780 \text{ cm}^{-1}$ which is not shown). **In the center of the spectrum is the Rayleigh scatter peak at the laser wavelength. This peak is millions of times more intense than the Raman scatter and is therefore normally blocked by a notch or edge filter in the Raman spectrometer but was included here for clarity.** Symmetrically placed on either side of the Rayleigh peak are the three Stokes and three Anti-Stokes peaks corresponding to the three most intense Raman active vibrations of CCl_4 . It can be seen that the Anti-Stokes lines are much weaker than the Stokes due to the larger population of molecules in the ground vibrational level of each mode. CCl_4 has one of the simplest Raman Spectra but the same principle applies for all samples: Raman spectroscopy is used to measure the unique vibrational fingerprint of the sample and from that information chemical, structural and physical properties can be determined






Raman spectra: What to do with them?

Raman spectroscopy is a powerful technique for analyzing the molecular structure and composition of materials. Raman spectra are obtained by irradiating a sample with a monochromatic laser beam and measuring the scattered light at different wavelengths. The scattered light contains information about the vibrational modes of the molecules in the sample, which can be used to identify and quantify them.

To interpret raman spectra, one needs to **compare the observed peaks with reference data from databases or literature sources**. **The peak positions, intensities, shapes, and widths can provide clues about the type, number, and arrangement of atoms and bonds in the sample**. Some common features of raman spectra are:

- The peak at zero wavenumber (also called Rayleigh scattering) corresponds to the elastic scattering of the laser light by the sample. It is usually much stronger than the other peaks and can be filtered out by using a notch filter or a monochromator.
- The peaks at positive and negative wavenumbers (also called Stokes and anti-Stokes scattering) correspond to the inelastic scattering of the laser light by the sample. They reflect the energy difference between the initial and final vibrational states of the molecules. The Stokes peaks are usually stronger than the anti-Stokes peaks because they involve a lower-energy transition.
- **The peak positions depend on the frequency of the vibrational modes of the molecules**, which are determined by their mass, bond strength, and symmetry. **Different functional groups have characteristic peak positions that can be used to identify them. For example, C-H stretching vibrations usually occur around 3000 cm⁻¹, C=O stretching vibrations around 1700 cm⁻¹, and C-C stretching vibrations around 1600 cm⁻¹.**
- The peak intensities depend on the number of molecules that contribute to a given vibrational mode, as well as their orientation and polarization relative to the laser beam. The intensity can also be affected by factors such as temperature, pressure, concentration, and sample preparation. Therefore, **peak intensities are not always reliable for quantitative analysis and may require calibration or normalization**.
- The **peak shapes and widths depend on the homogeneity and complexity of the sample**, as well as the resolution and sensitivity of the instrument. Broad peaks indicate a large distribution of vibrational frequencies, which can be caused by disorder, defects, or interactions among molecules. Narrow peaks indicate a well-defined vibrational frequency, which can be caused by crystallinity, purity, or isolation of molecules.

Raman spectra




Uses of Raman spectroscopy

Raman spectroscopy is a powerful analytical technique that can provide information about the molecular structure, composition, and interactions of various materials.

Raman spectroscopy uses the scattering of light by molecules to detect their vibrational modes, which are characteristic of their chemical bonds and symmetry.

Raman spectroscopy has many applications in different fields, such as chemistry, physics, biology, medicine, geology, and engineering. Some of the uses of Raman spectroscopy are:

- **Identification and characterization** of organic and inorganic compounds, such as drugs, polymers, minerals, and crystals.
- **Analysis of molecular interactions and conformations**, such as hydrogen bonding, solvation, and folding.
- **Detection and quantification of impurities, contaminants, and adulterants** in samples, such as food, pharmaceuticals, and environmental samples.
- Investigation of **phase transitions, crystallization, and polymorphism** in materials, such as metals, alloys, ceramics, and polymers.
- Study of **surface phenomena and nanostructures**, such as catalysis, corrosion, adsorption, and plasmonics.
- **Imaging and mapping of biological tissues** and cells, such as skin, bone, cartilage, blood vessels, and neurons.



Infrared vs Raman spectroscopy

Some of the **advantages** of Raman spectroscopy over infrared spectroscopy are:

- Raman spectroscopy can be used with **any wavelength** of light, whereas infrared spectroscopy requires a specific range of wavelengths that correspond to the vibrational frequencies of the molecule.
- Raman spectroscopy can be performed in **aqueous solutions**, whereas infrared spectroscopy is hindered by the strong absorption of water.
- Raman spectroscopy can analyze samples with **minimal or no sample preparation**, whereas infrared spectroscopy may require grinding, pressing, or diluting the sample.
- Raman spectroscopy can provide **information about crystal structure and orientation**, whereas infrared spectroscopy cannot.


Some of the **disadvantages** of Raman spectroscopy compared to infrared spectroscopy are:

- Raman spectroscopy has a **lower signal-to-noise ratio** than infrared spectroscopy, which means it requires more powerful lasers and longer acquisition times.
- Raman spectroscopy may **cause damage to the sample** due to heating or photodegradation from the laser beam, whereas infrared spectroscopy is less invasive.
- Raman spectroscopy **may suffer from fluorescence interference** from some organic compounds or impurities, whereas infrared spectroscopy does not.
- Raman spectroscopy may require calibration and normalization procedures to **account for variations in laser power** and detector sensitivity, whereas infrared spectroscopy does not.

However, they complement each other

Refractometry


Refractometry is a technique used to measure the refractive index of a substance. **The refractive index is a measure of how much a substance bends light as it passes through it.** By measuring the angle at which light is bent as it passes through a sample, the refractive index can be determined. Refractometry is commonly used in chemistry, biology, and other fields to determine the composition and purity of substances. It is also used in the food industry to measure the sugar content of fruits and other products.



Differences between refractometers and spectrophotometers


A **refractometer** measures the bending of light as it passes through a liquid or transparent solid, while a **spectrophotometer** measures the amount of light absorbed or transmitted by a sample at different wavelengths.

Refractometers are commonly used to measure the concentration of a solution based on its refractive index, while spectrophotometers are used to determine the chemical composition of a sample by analyzing its absorption or transmission spectrum.



Refractive index (RI)

The refractive index (RI) is a fundamental optical property that describes the way light propagates through a medium. It is defined as the ratio of the speed of light in a vacuum to the speed of light in a given material. **The RI of a material depends on its chemical composition and physical properties such as temperature, pressure, and density.** In general, the higher the RI of a material, the more it bends light when it passes through it. The refractive index is an important parameter in a wide range of scientific and industrial applications, including optics, chemistry, materials science, and food science.



Refractometry

The refractive index (RI) of a substance is affected by several parameters, including temperature, pressure, concentration, and wavelength of light.

Temperature: As temperature increases, the RI of a substance generally decreases. This is because the speed of light in a medium increases with temperature, which reduces the bending of light and hence the RI.

Pressure: Pressure can also affect the RI of a substance, but the effect is generally small. Higher pressures can cause slight changes in the density of a substance, which can alter the RI.

Concentration: The RI of a substance is directly proportional to its concentration. This property is used in refractometry to determine the concentration of a solute in a solution.

Wavelength of light: The RI of a substance varies with the wavelength of light used. This property is exploited in spectral refractometry, which uses different wavelengths of light to determine the RI of a substance at each wavelength.


The modification of the uses of refractometry depends on the parameter being varied. For example, changes in temperature or pressure may affect the accuracy of RI measurements, **while changes in concentration can provide a means for determining the concentration of a solute in a solution.** Varying the wavelength of light used can provide more detailed information about the chemical composition of a substance.

Brix degree

YIERYI

Brix% : 28-62

A Brix degree is a unit of measurement used to determine the sugar content in a liquid, typically in fruit juices, wine, and other beverages. It is named after Adolf Brix, who developed the method in the late 19th century. **One degree Brix (°Bx) is equal to 1 gram of sucrose in 100 grams of solution.** For example, a solution with a Brix value of 10 contains 10% sugar by weight at 20°C. The Brix degree is used in the food and beverage industry to determine the quality and ripeness of fruits and vegetables, as well as to monitor the fermentation process in winemaking and brewing.



Refractometry: Uses (Examples)

Measuring the concentration of a solution: Refractometers are commonly used in the food and beverage industry to measure the sugar content of fruits, vegetables, and juices. They can also be used to measure the concentration of other substances such as salt, oil, and alcohol.

Testing the purity of a substance: Refractometers can be used to test the purity of a substance by measuring its refractive index. This is particularly useful in the pharmaceutical industry where it is important to ensure that drugs are pure and free from contamination.

Checking the quality of lubricants and oils: Refractometers can be used to check the quality of lubricants and oils by measuring their viscosity and density. This is important in the automotive and aviation industries where the performance of engines depends on the quality of the lubricants and oils used.


Monitoring the concentration of coolant in a car: Refractometers can be used to monitor the concentration of coolant in a car. This is important as the coolant helps to regulate the temperature of the engine and prevent it from overheating.

Testing the salinity of water: Refractometers can be used to test the salinity of water in aquariums and marine environments. This is important as different species of fish and marine life require different levels of salinity to survive.





Differential Scanner Calorimetry (DSC)

DSC (differential scanning calorimetry) is an analytical technique that measures the **thermal behavior of a sample as a function of temperature or time**. It is commonly used in the field of materials science to study the thermal properties of polymers, metals, ceramics, and other materials. In DSC, a sample and a reference material are placed in separate cells in a calorimeter. The cells are heated or cooled at a controlled rate, and the heat flow between the cells is measured. As the sample undergoes a thermal transition, such as melting or crystallization, it absorbs or releases heat, causing a change in the heat flow signal. This signal is plotted as a function of temperature or time, generating a DSC curve that provides information about the thermal properties of the sample.



DSC: Principle

The principle of differential scanning calorimetry (DSC) is based on the **measurement of the difference in heat flow between a sample and a reference material as they are subjected to a controlled temperature program**. As the temperature of the sample is changed, it undergoes thermal transitions such as melting, crystallization, or glass transitions, which are accompanied by changes in heat capacity. These changes in heat capacity result in a difference in the heat flow between the sample and the reference material, which is detected by a thermocouple and recorded as a function of temperature. By analyzing the resulting DSC curve, information about the thermal properties of the sample, such as its melting point, glass transition temperature, and heat of fusion, can be obtained.



DSC: Uses

Characterization of thermal properties: DSC is used to determine the thermal properties of materials such as melting point, glass transition temperature, crystallization temperature, and heat capacity.

Quality control: DSC is used in quality control to ensure that materials meet specific standards and to detect any impurities or defects.

Material identification: DSC can be used to identify unknown materials by comparing their thermal properties to those of known materials.

Stability studies: DSC is used to study the stability of materials under different conditions such as temperature, humidity, and light exposure.

Formulation development: DSC is used in the development of new formulations for pharmaceuticals, cosmetics, and other products to optimize their thermal properties.

Polymer analysis: DSC is commonly used in the analysis of polymers to determine their thermal properties, such as glass transition temperature and melting point, which are important for understanding their behavior in different applications.

